Mathis
AIR BOY
Becherfärbeapparat
Typ AB

Mathis
Mathis
Mathis
Mathis
Mathis

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General

The AIR BOY is a beaker dyeing unit with a heating system without liquid heating medium. Already more than 10 years ago Mathis AG has introduced the well known IR heating system in larger numbers with the beaker dyer LABOMAT. This unit is perfectly suitable for dyeings under HT conditions. The AIR BOY is now the corresponding link for dyeings in the atmospheric temperature range with 2 independent sets of dye beakers with 6 or 12 beakers each. For the heating of the unit hot air is used as a heat source. The air circulation is made via an efficient fan. The temperature cycle for both set of dye beakers is controlled via the same UNIPROGRAMMER whereby both cycles can run completely independent of each other. In order to have the highest precision during the dye cycle the bath temperature is measured directly in one of the beakers per set. Visual access during the dye cycle is enabled by the full view bath container made of Borsilicate glass.

The AIR BOY needs a space of approx. 0.8 m² and is available in two versions for 2 x 6 or 2 x 12 beakers.

Since the AIR BOY is a beaker dyer for the atmospheric temperature range the dye beakers are always easily accessible. Due to this reason the apparatus is specially suitable for reactiv dyeings or other dye processes where chemicals have to be added. This can be done by means of a hand syringe or a multiple dosing pump. With a simple modification set the use of beakers with a height of 240 mm or 400 mm is possible.

Sample movement

The skeins, piece goods or loose material are being moved up and down in the dye liquor by means of a synchronized stroke movement.

The stroke frequency, within 15–60 stroke movements per minute, can be entered at the UNIPROGRAMMER. The stroke length is 40 mm.

Loading/unloading

For loading and unloading the flat-top section including material holders can be raised up or lowered electrically. This enables an easy loading of the samples and unloading after the dyeing process. At the end of the dyeing time the beakers are cooled by means of air circulation. The cooling efficiency is increased by a water cooling coil. This principle avoids an outlet of hot air into the laboratory and does not influence the laboratory environment.

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Program control


The programmable process controller Mathis UNIPROGRAMMER has a memory capacity for 70 user programs, 30 steps each. The process controller is very user-friendly and simple to use, because a lot of information is displayed in text form. Beside automatic program execution, also a simple temperature can be entered with its corresponding gradient plus following holding time.

The programs are entered by using dialog mode and menu technique with english text display. Various function outputs as well as a selection of system commands enable a large variety of individual program formation.

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On a clearly visible display there are two possibilities of indication, either the dye curve of the running program with current position indicator or the corresponding program details in text mode. The actual value for temperature as well as status indication for heating or cooling for both set of dye beakers are constantly displayed.

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Programming is done either via keyboard in dialog mode with the process controller or it can be done on a personal computer with corresponding software. All program data, machine parameters as well as further information are memorized on a memory card. By the use of this memory card the data can be transferred to other controllers of the same type. The same controller enables also a great number of interlockings of various machine functions. This reduces operating errors to a great extent.

Also a great number of control functions can be realised for example temperature sensor control etc. which enables a perfect reliability. In case of errors which might happen during the program execution or in case of machine errors an alarm sounds and/or the alarm unit type FB is intermittent flashing if it is connected as well as a corresponding indication in text form will be displayed.

Protocol printer

A thermal printer with 40 digits per line is available for printing of a protocol. The print out can be used as a tool for quality control of the actually processed dye program including any errors which might have occurred.

Process control via personal computer (PC)

If this option is required the data can be transferred to the personal computer by the memory card. It is also possible to connect the controller directly to the PC in order to create or to alter the programs at a central computer as well as to get actual program information of the running program. In this combination also a print-out of the finished program (e.g. for ISO 9001) can be made in which also program stops and possible error statements are listed.

Dosing

There are three possibilities for the dosing:

1. By means of a micro pipette

With this method the required quantity is added into the dye beaker by means of a micro pipette. After the adding the same quantity for the next beaker is automatically prepared. Then dosing is done in the following beaker. Following sizes of micro pipettes are available: 0–5 or 0–10 ml. Larger quantities can, where required, also be obtained by multiple injection.

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2. Dosing by means of programmable diluter

This process is very similar to the one mentioned above, however the required dosing quantity can be programmed. After inserting the needle into the beaker the adding is done by pressing a button. Then introduce the needle into the next beaker and dose again. The diluter is equipped with a standard injector of 0–10 ml. Larger quantities can be obtained by multiple injection. For dosing of several products, a multi syringe diluter is available on request.

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3. Multichannel dosing pump

For the continuous, programmable dosing in modern processes with reactive dyes. With suitable software following dosing curves are possible progressive, linear and degressive whereby dosing quantity and time can be programmed.

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Werner Mathis AG
Werner Mathis U.S.A. Inc.
Concord N.C. 28026, U.S.A.